

Product Datasheet

Human Immunoglobulin M (IgM) Antibody HM098

Product Name	Human Immunoglobulin M (IgM) Antibody
Catalog Number	HM098
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Human IgM
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Human IgM
Preservatives	0.1% Proclin 300