

## Product Datasheet

### Estrone-3-Sulfate (E1S) Antibody HM1016

<b>Product Name</b>	Estrone-3-Sulfate (E1S) Antibody
<b>Catalog Number</b>	HM1016
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Estrone-3-Sulfate-KLH
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90% by SDS-PAGE
<b>Buffer</b>	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
<b>Cross Reactivity</b>	Estriol <6%; Estriol-3-glucuronide, 17 $\beta$ -Estradiol, Estradiol-3-glucuronide <5%; 17 $\alpha$ -Ethinylestradiol <3%; Estrone-3-glucuronide, 17 $\alpha$ -hydroxyprogesterone, Progesterone, Pregnenediol-3-glucuronide, Testosterone, Aldosterone, Androstenedione <1%
<b>Specificity</b>	Estrone and Estrone-3-sulfate
<b>Matched Pairs</b>	LA550 or LA551 (capture)
<b>Preservatives</b>	0.05% Sodium Azide