

Product Datasheet

Insulin-Like Growth Factor Binding Protein 1 (IGFBP-1) Antibody HM1024

Product Name	Insulin-Like Growth Factor Binding Protein 1 (IGFBP-1) Antibody
Catalog Number	HM1024
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Human IGFBP-1
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	IGFBP-1
Matched Pairs	HM1025 (detection)
Preservatives	0.09% Sodium Azide