

Product Datasheet

African Swine Fever Virus (ASFV-p30) Antibody HM1027

Product Name	African Swine Fever Virus (ASFV-p30) Antibody
Catalog Number	HM1027
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant ASFV strain Morara/Georgia
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1031 (capture)
Preservatives	0.05% Sodium Azide