

Product Datasheet

African Swine Fever Virus (ASFV-p30) Antibody HM1029

| | |
|----------------------------|---|
| Product Name | African Swine Fever Virus (ASFV-p30) Antibody |
| Catalog Number | HM1029 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant ASFV strain Morara/Georgia |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Matched Pairs | HM1031 (capture) |
| Preservatives | 0.05% Sodium Azide |