

Product Datasheet

Pancreatic Elastase Antibody HM1033

Product Name	Pancreatic Elastase Antibody
Catalog Number	HM1033
Species/Host	Mouse
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with pig elastase; pig, bovine, human hemaglobin; pig, bovine, human transferrin; human calprotectin; bovine, human lactoferrin
Specificity	Pancreatic Elastase
Matched Pairs	HM1034 (capture)
Preservatives	0.09% Sodium Azide