

Product Datasheet

Hydrocodone Antibody HM1037

Product Name	Hydrocodone Antibody
Catalog Number	HM1037
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Hydrocodone-BSA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Detection limit = 10 ng/ml.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Codeine, Oxycodone 30%; Hydromorphone, Oxymorphone <2%; 6-MAM, Morphine <1%
Specificity	Hydrocodone
Matched Pairs	LA289 (capture)
Preservatives	0.05% Sodium Azide