

Product Datasheet

Pancreatic Lipase Antibody, Canine HM1048

Product Name	Pancreatic Lipase Antibody, Canine
Catalog Number	HM1048
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Native Canine Pancreatic Lipase
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Matched Pairs	HM1049 (capture or detection)
Preservatives	0.02% Sodium Azide