

## Product Datasheet

### Coronavirus (COVID-19 & SARS-CoV NP) Antibody HM1054

<b>Product Name</b>	Coronavirus (COVID-19 & SARS-CoV NP) Antibody
<b>Catalog Number</b>	HM1054
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, WB, IFA
<b>Immunogen</b>	Recombinant SARS-CoV (aa 1-49)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	No reactivity with Canine Coronavirus, Feline Coronavirus (FIP-1 & 2), TGEV, Coronavirus OC43 & 229E
<b>Specificity</b>	SARS-CoV-2 (COVID-19) & SARS-CoV Nucleoprotein (NP)
<b>Preservatives</b>	0.1% Sodium Azide