

Product Datasheet

Coronavirus (COVID-19, MERS & SARS-CoV NP) Antibody HM1058

Product Name	Coronavirus (COVID-19, MERS & SARS-CoV NP) Antibody
Catalog Number	HM1058
Species/Host	Mouse
Tested Applications	ELISA, IFA
Immunogen	Recombinant Coronavirus nucleoprotein peptide
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with Coronavirus OC43, NL63, 229E
Specificity	SARS-CoV-2 (COVID-19), MERS & SARS-CoV Nucleoprotein (NP)
Preservatives	0.1% Sodium Azide