

Product Datasheet

Hepatitis Be Antigen (HBeAg) Antibody HM1061

Product Name	Hepatitis Be Antigen (HBeAg) Antibody
Catalog Number	HM1061
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	LA577 (capture)