

## Product Datasheet

### Coronavirus (COVID-19 NP) Antibody HM1069

<b>Product Name</b>	Coronavirus (COVID-19 NP) Antibody
<b>Catalog Number</b>	HM1069
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Native Coronavirus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Use with native SARS-CoV-2 for best results
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95%
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	Canine Coronavirus
<b>Specificity</b>	SARS-CoV-2 (COVID-19) & SARS-CoV Nucleoprotein (NP)
<b>Matched Pairs</b>	HM1068
<b>Preservatives</b>	0.09% Sodium Azide