

## Product Datasheet

### Human Immunoglobulin G (IgG) Antibody HM1070

<b>Product Name</b>	Human Immunoglobulin G (IgG) Antibody
<b>Catalog Number</b>	HM1070
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Human IgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with Human IgA, IgM, IgE
<b>Specificity</b>	Human IgG (IgG1, IgG2a, IgG3, IgG4)
<b>Preservatives</b>	0.09% Sodium Azide