

Product Datasheet

Human Immunoglobulin G (IgG) Antibody HM1070

Product Name	Human Immunoglobulin G (IgG) Antibody
Catalog Number	HM1070
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Human IgG
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Human IgA, IgM, IgE
Specificity	Human IgG (IgG1, IgG2a, IgG3, IgG4)
Preservatives	0.09% Sodium Azide