

## Product Datasheet

### Coronavirus (COVID-19 NP) Antibody HM1074

<b>Product Name</b>	Coronavirus (COVID-19 NP) Antibody
<b>Catalog Number</b>	HM1074
<b>Species/Host</b>	Cultured cells
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	Humanized IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95%
<b>Buffer</b>	20 mM Tris with 50 mM glycine
<b>Matched Pairs</b>	HM1075