

Product Datasheet

Coronavirus (COVID-19 NP) Antibody HM1075

Product Name	Coronavirus (COVID-19 NP) Antibody
Catalog Number	HM1075
Species/Host	Cultured cells
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	Humanized IgG
Clonality	Monoclonal
Purity	>95%
Buffer	Phosphate Buffered Saline
Matched Pairs	HM1074