

Product Datasheet

Coronavirus (COVID-19 NP) Antibody HM1077

Product Name	Coronavirus (COVID-19 NP) Antibody
Catalog Number	HM1077
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	SARS-CoV-2 (COVID-19) nucleoprotein
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	Carbonate-Bicarbonate Buffer, pH 7.5
Cross Reactivity	SARS-CoV
Specificity	SARS-CoV-2 (COVID-19)