

## Product Datasheet

### Human CD62E ELISA Kit 850.530

<b>Product Name</b>	Human CD62E ELISA Kit
<b>Catalog Number</b>	850.530
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	1-32 ng/ml
<b>Sensitivity</b>	0.5 ng/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	E-selectin / ELAM-1
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human adhesion molecules
<b>Specificity</b>	Recognizes both natural and recombinant human soluble ELAM-1
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant