

Product Datasheet

Synuclein, alpha (?-Syn) Antibody HM1089

Product Name	Synuclein, alpha (α -Syn) Antibody
Catalog Number	HM1089
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant human α -Synuclein (insect cell expressed)
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: \leq -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1093 or HM1094 (detection)
Preservatives	0.05% Sodium Azide