

Product Datasheet

Neurofilament L (NfL) Antibody HM1095

Product Name	Neurofilament L (NfL) Antibody
Catalog Number	HM1095
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant NfL protein (insect cell expressed)
Alternative Names	NfL Neurofilament Light Chain Neurofilament Light Peptide
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Pair with HM1096 or HM1184 for best results in Lateral Flow
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1096, HM1098, HM1101 or HM1184 (detection)
Preservatives	0.05% Sodium Azide