

Product Datasheet

Neurofilament L (NfL) Antibody HM1095

| | |
|----------------------------|---|
| Product Name | Neurofilament L (NfL) Antibody |
| Catalog Number | HM1095 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant NfL protein (insect cell expressed) |
| Alternative Names | NfL Neurofilament Light Chain Neurofilament Light Peptide |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Note | Pair with HM1096 or HM1184 for best results in Lateral Flow |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Matched Pairs | HM1096, HM1098, HM1101 or HM1184 (detection) |
| Preservatives | 0.05% Sodium Azide |