

Product Datasheet

Neurofilament L (NfL) Antibody HM1097

Product Name	Neurofilament L (NfL) Antibody
Catalog Number	HM1097
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Recombinant NfL protein (insect cell expressed)
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1096 or HM1101 (capture)
Preservatives	0.05% Sodium Azide