

Product Datasheet

Neurofilament L (NfL) Antibody HM1098

| | |
|----------------------------|---|
| Product Name | Neurofilament L (NfL) Antibody |
| Catalog Number | HM1098 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant NfL protein (insect cell expressed) |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥90% |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Matched Pairs | HM1095 (capture) |
| Preservatives | 0.05% Sodium Azide |