

Product Datasheet

Neurofilament L (NfL) Antibody HM1099

Product Name	Neurofilament L (NfL) Antibody
Catalog Number	HM1099
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant NfL protein (insect cell expressed)
Alternative Names	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS), Amyotrophic lateral sclerosis (ALS)
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Preservatives	0.05% Sodium Azide