

## Product Datasheet

### Neurofilament L (NfL) Antibody HM1100

<b>Product Name</b>	Neurofilament L (NfL) Antibody
<b>Catalog Number</b>	HM1100
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant NfL protein (insect cell expressed)
<b>Alternative Names</b>	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS), Amyotrophic lateral sclerosis (ALS)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
<b>Matched Pairs</b>	HM1096 or HM1101 (capture)
<b>Preservatives</b>	0.05% Sodium Azide