

Product Datasheet

Coronavirus (SARS-CoV NP) Antibody HM1127

Product Name	Coronavirus (SARS-CoV NP) Antibody
Catalog Number	HM1127
Species/Host	Mouse
Tested Applications	LF
Immunogen	Recombinant SARS-CoV-2 (COVID-19) nucleoprotein
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with SARS-CoV-2 (COVID-19)
Specificity	SARS-CoV
Matched Pairs	HM1126 (detection)
Preservatives	0.05% Sodium Azide