

Product Datasheet

Cancer Antigen 15-3 (CA 15-3) MUC 1 Antibody HM1131

Product Name	Cancer Antigen 15-3 (CA 15-3) MUC 1 Antibody
Catalog Number	HM1131
Species/Host	CHO Cells
Tested Applications	CLIA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffer, pH 7.4 with 0.15M NaCl
Matched Pairs	HM1130 (detection)
Preservatives	0.1% Proclin 300