

Product Datasheet

Fentanyl Antibody HM1133

Product Name Fentanyl Antibody

Catalog Number HM1133

Species/Host Mouse

Tested Applications LF, ELISA

Immunogen Fentanyl-KLH

Brand EastCoast Bio

Form/Appearance Ion Exchange Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1

Clonality Monoclonal

Purity >90%

Buffer 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride

Cross Reactivity Acetyl Norfentanyl 90%; Acetyl Fentanyl & Norfentanyl 60%; Alfentanil 30%; Carfentanil 1%

Specificity Fentanyl

Matched Pairs LA395 (capture)

Preservatives 0.05% Sodium Azide