

Product Datasheet

Pepsin Antibody HM1136

Product Name	Pepsin Antibody
Catalog Number	HM1136
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Human Pepsin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Matched Pairs	HM1134 (detection)
Preservatives	0.09% Sodium Azide