

Product Datasheet

Estrone-3-Glucuronide (E1G) Antibody HM1150

Product Name	Estrone-3-Glucuronide (E1G) Antibody
Catalog Number	HM1150
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Estrone-3-Glucuronide-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Lateral Flow linear range 0.5-10 ng/ml
Isotype	IgG1
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Cross Reactivity	Estrone 60%; Estriol <6%; Estradiol-3-glucuronide, 17β-Estradiol, Estriol-3-Glucuronide <5%; 17α-Ethinylestradiol <3%; 17α-hydroxyprogesterone, Testosterone, Progesterone, Pregnanediol-3-Glucuronide, Aldosterone, Androstenedione <1%
Specificity	Estrone-3-glucuronide
Matched Pairs	LA655 or LA669 (capture)
Preservatives	0.05% Sodium Azide