

## Product Datasheet

### Risperidone Antibody HM1154

<b>Product Name</b>	Risperidone Antibody
<b>Catalog Number</b>	HM1154
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Risperidone-BSA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90% by SDS-PAGE
<b>Buffer</b>	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	Amitriptyline, Aripiprazole, Clomipramine, Clonidine, Cyclobenzaprine, Desipramine, Doxepin, Imipramine, Maprotiline, Nordoxepin, Nortriptyline, Olanzapine, Quetiapine, Trimipramine, Ziprasidone <1%
<b>Specificity</b>	Risperidone
<b>Matched Pairs</b>	LA656 (capture)
<b>Preservatives</b>	0.05% Sodium Azide