

Product Datasheet

Mouse Immunoglobulin (IgG) Antibody HM1161

Product Name	Mouse Immunoglobulin (IgG) Antibody
Catalog Number	HM1161
Species/Host	Mouse
Tested Applications	LF, Blocking
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl
Preservatives	0.1% Proclin 300