

Product Datasheet

Influenza B Virus Antibody HM1165

Product Name	Influenza B Virus Antibody
Catalog Number	HM1165
Species/Host	CHO Cells
Tested Applications	LF
Immunogen	Inactivated Influenza B Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza A, Adenovirus, RSV, Parainfluenza, M. pneumoniae or SARS-CoV-2
Specificity	Influenza B Virus Nucleoprotein
Matched Pairs	HM1166 (detection)
Preservatives	0.1% Proclin 300