

Product Datasheet

Coronavirus (COVID-19 NP) Antibody, Rabbit Recombinant HM1180

Product Name	Coronavirus (COVID-19 NP) Antibody, Rabbit Recombinant
Catalog Number	HM1180
Species/Host	Rabbit
Tested Applications	LF
Immunogen	Recombinant SARS-CoV-2 (COVID-19) nucleoprotein
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Recombinant Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	High binding affinity and reactivity with native and recombinant antigens. Highest sensitivity when used in combination.
Isotype	IgG
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with MERS, Coronavirus NL63, 229E
Specificity	SARS-CoV-2 (COVID-19) & SARS-CoV Nucleoprotein (NP)
Matched Pairs	HM1181 or HM1183 (detection)
Preservatives	0.05% Sodium Azide