

## Product Datasheet

### Human Chorionic Gonadotropin $\beta$ (HCG $\beta$ ) Antibody HM119

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Human Chorionic Gonadotropin $\beta$ (HCG $\beta$ ) Antibody |
| <b>Catalog Number</b>      | HM119  |
| <b>Species/Host</b>        | Mouse  |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Human Chorionic Gonadotropin (beta subunit)                  |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody                                 |
| <b>Storage</b>             | 2-8°C  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Purity</b>              | >95%   |
| <b>Buffer</b>              | 0.05 M Phosphate Buffered Saline, pH 7.4                     |
| <b>Cross Reactivity</b>    | $\beta$ hCG 100%, ihCG 115%, LH < 1%, FSH < 1%, TSH < 1%     |
| <b>Specificity</b>         | HCG $\beta$ subunit  |
| <b>Matched Pairs</b>       | HM118 (capture)  |
| <b>Preservatives</b>       | 0.1% Sodium Azide  |