

Product Datasheet

Prolactin (PRL) Antibody HM120

Product Name	Prolactin (PRL) Antibody
Catalog Number	HM120
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Human Prolactin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG
Clonality	Monoclonal
Buffer	0.05M Phosphate Buffered Saline, pH 7.4
Specificity	Human Prolactin
Matched Pairs	FA310 or FA312 (detection)
Preservatives	0.1% Sodium Azide