

Product Datasheet

Lysergic acid diethylamide (LSD) Antibody HM1201

Product Name	Lysergic acid diethylamide (LSD) Antibody
Catalog Number	HM1201
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	LSD-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	PBS, pH 7.2
Specificity	LSD
Matched Pairs	LA783 (capture)
Preservatives	0.05% Sodium Azide