

Product Datasheet

Barbiturate Antibody HM137

Product Name Barbiturate Antibody

Catalog Number HM137

Species/Host Mouse

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1, kappa

Clonality Monoclonal

Purity $\geq 98\%$ by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Preservatives 0.05% Sodium Azide