

## Product Datasheet

### Morphine Antibody HM139

**Product Name** Morphine Antibody

**Catalog Number** HM139

**Species/Host** Mouse

**Tested Applications** LF, ELISA, IFA

**Brand** EastCoast Bio

**Form/Appearance** Protein A Purified Monoclonal Antibody

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Isotype** IgG1, lamda

**Clonality** Monoclonal

**Purity** >98% by SDS-PAGE

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.2

**Specificity** Morphine

**Preservatives** 0.05% Sodium Azide