

## Product Datasheet

### Gentamicin Antibody HM141

**Product Name** Gentamicin Antibody

**Catalog Number** HM141

**Species/Host** Mouse

**Tested Applications** LF, ELISA, IFA

**Brand** EastCoast Bio

**Form/Appearance** Protein A Purified Monoclonal Antibody

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Isotype** IgG1, kappa

**Clonality** Monoclonal

**Purity** >98%

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.2

**Preservatives** 0.05% Sodium Azide