

Product Datasheet

Phenytoin (Dilantin) Antibody HM142

Product Name Phenytoin (Dilantin) Antibody

Catalog Number HM142

Species/Host Mouse

Tested Applications LF, ELISA, IFA

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1, kappa

Clonality Monoclonal

Purity >98%

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Preservatives 0.05% Sodium Azide