

Product Datasheet

Beta-2 Microglobulin (B2M) Antibody HM143

Product Name	Beta-2 Microglobulin (B2M) Antibody
Catalog Number	HM143
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>98%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservatives	0.05% Sodium Azide