

Product Datasheet

NT-proBNP (38-44) Antibody HM146

Product Name	NT-proBNP (38-44) Antibody
Catalog Number	HM146
Species/Host	Mouse
Tested Applications	LF, CLIA
Immunogen	Native Human ProBNP antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10mM Phosphate Buffered Saline, pH 7.4
Specificity	NT-proBNP
Matched Pairs	HM148 (detection)
Preservatives	0.1% Proclin 300