

## Product Datasheet

### Methamphetamine Antibody HM165

<b>Product Name</b>	Methamphetamine Antibody
<b>Catalog Number</b>	HM165
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	No cross reactivity with Amphetamine
<b>Specificity</b>	Methamphetamine
<b>Matched Pairs</b>	LA201 (capture)
<b>Preservatives</b>	0.02% Sodium Azide