

Product Datasheet

Parainfluenza Virus type 3 (PIV-3) Antibody HM181

Product Name	Parainfluenza Virus type 3 (PIV-3) Antibody
Catalog Number	HM181
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Bovine PIV-3
Specificity	Parainfluenza Virus type 3 hemagglutinin
Matched Pairs	HM181 (self-pairs) or HM182
Preservatives	0.1% Sodium Azide