

Product Datasheet

Streptococcus pneumoniae Antibody HM187

Product Name	Streptococcus pneumoniae Antibody
Catalog Number	HM187
Species/Host	Rabbit
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	S. pneumoniae common cell wall polysaccharide (urinary antigen)
Matched Pairs	HM188 (capture)
Preservatives	0.1% Sodium Azide