

Product Datasheet

Troponin I (41-49) Antibody HM192

| | |
|----------------------------|--|
| Product Name | Troponin I (41-49) Antibody |
| Catalog Number | HM192 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA, WB |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | 2-8°C |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Cross Reactivity | Does not cross-react with skeletal muscle Troponin I |
| Specificity | Troponin I (epitope 41-49 a.a.) |
| Preservatives | 0.09% Sodium Azide |