

Product Datasheet

Troponin I (13-22) Antibody HM200

Product Name	Troponin I (13-22) Antibody
Catalog Number	HM200
Species/Host	Mouse
Tested Applications	LF, WB
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG3
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	Does not cross-react with skeletal muscle Troponin I
Specificity	Troponin I (epitope 13-22 a.a.)
Matched Pairs	HM197 (detection)
Preservatives	0.09% Sodium Azide