

Product Datasheet

NT-proBNP (39-76) Antibody HM209

Product Name NT-proBNP (39-76) Antibody

Catalog Number HM209

Species/Host Mouse

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Purified Monoclonal Antibody

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1

Clonality Monoclonal

Purity >90%

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Specificity NT-proBNP

Matched Pairs HM210 (detection)

Preservatives 0.1% Proclin 300