

Product Datasheet

NT-proBNP (13-27) Antibody HM210

Product Name	NT-proBNP (13-27) Antibody
Catalog Number	HM210
Species/Host	Mouse
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥90%
Buffer	10 mM Phosphate Buffer, pH 7.4 with 0.15M NaCl
Specificity	NT-proBNP
Matched Pairs	HM209 or HM899 (capture)
Preservatives	0.1% Proclin 300